

# POLYGLYCAN®

**(Hyaluronic Acid + Chondroitin sulfates C4 & C6 + N-acetyl-D-glucosamine Sterile Solution)**

**STERILE POST-SURGICAL LAVAGE, SYNOVIAL FLUID REPLACEMENT, VISCOSUPPLEMENT**

**Veterinary Use Only  
FOR USE IN HORSES**

Patent Nos.: 6,979,679 / 7,485,629 / 8,455,458



## INDICATIONS

Polyglycan® is indicated for the replacement of lost or damaged synovial fluid. Polyglycan® contains naturally occurring components of the synovia that play a central part in maintaining the homeostatic environment of the joint.

## BENEFITS

- **Helps replenish** lost synovial fluid with naturally occurring elements
- **Assists** in the normalization of synovial fluid viscosity
- **Post surgical lavage** of joint cavity and surfaces
- **Patented formulation**

## CAUTION

This device is restricted to use by or on the order of a licensed veterinarian.

## PACKAGING

LIST NO.	UNIT PACKAGE	CASE SIZE
1POL007	6 x 10 mL vials	40

**Bimeda, Inc.**  
One Tower Lane, Suite 2250  
Oakbrook Terrace, IL 60181 USA  
Toll Free Tel. 888-524-6332 • Toll Free Fax. 877-888-7035  
Email: [Info@BimedaUS.com](mailto:Info@BimedaUS.com)

See reverse side for additional information



# POLYGLYCAN®

*(Hyaluronic Acid + Chondroitin sulfates C4 & C6 + N-acetyl-D-glucosamine Sterile Solution)*

**STERILE POST-SURGICAL LAVAGE, SYNOVIAL FLUID REPLACEMENT, VISCOSUPPLEMENT**

***Veterinary Use Only  
FOR USE IN HORSES***

Patent Nos.: 6,979,679 / 7,485,629 / 8,455,458

**CAUTION:**

This device is restricted to use by or on the order of a licensed veterinarian.

**DESCRIPTION:**

Polyglycan® is a patented formulation designed to replace lost or damaged synovial fluid. Polyglycan® contains naturally occurring components of the synovia that play a central role in maintaining the homeostatic environment of the joint.

**COMPOSITION:**

Highly viscous aqueous solution of defined fractions of purified Hyaluronic acid, Chondroitin sulfates C4 & C6 in a 10% solution of N-acetyl-D-glucosamine.

**EACH VIAL CONTAINS:**

Hyaluronic acid sodium salt.....	50 mg
Sodium chondroitin sulfate.....	1000 mg
N-acetyl-D-glucosamine.....	1000 mg

**STORAGE:**

Refrigerate whenever possible. DO NOT FREEZE. Maintain product at or below 65°F. (18°C.) during shipment and storage.

**PROPERTIES:**

Glycosaminoglycans are important components of all extracellular tissue structures including cartilage and synovial fluid. The active components in Polyglycan® exhibit viscoelastic and polyionic properties similar in nature to synovial fluid.

Introduce Polyglycan® into the synovial space as a post-surgical lavage, synovial fluid replacement or viscosupplement to assist in the normalization of the joint.

**DIRECTIONS:**

Instill Polyglycan® into the synovial space to replenish lost synovial fluid. Quantity to be determined by Veterinarian Surgeon. Discard any unused portion of vial contents.

**WARNING:**

**Do not administer to animals that are to be slaughtered for use in food. Keep out of reach of children.**